

Product Datasheet

Human Sclerostin Domain-containing Protein 1 ELISA Kit (orb1210058)

Description	Human Sclerostin Domain-containing Protein 1, SOSTDC1 ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	SOSTDC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,307
Uniprot ID	Q6X4U4
Solubility (25°C)	48ng/ml
Sensitivity	0.06ng/ml
Entrez	25928
Expiration Date	Please enquire.