

Product Datasheet

Human Mesencephalic Astrocyte-derived Neurotrophic Factor ELISA Kit (orb1210062)

Description	Human Mesencephalic Astrocyte-derived Neurotrophic Factor, MANF ELISA Kit
Reactivity	Human
Range	10-3800ng/L
Tested Applications	ELISA
Target	MANF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,700
Uniprot ID	P55145
Solubility (25°C)	4000ng/L
Sensitivity	5.39ng/L
Entrez	7873
Expiration Date	Please enquire.