

Product Datasheet

Human Stress-70 Protein ELISA Kit (orb1210064)

Description	Human Stress-70 Protein, Mitochondrial, HSPA9 ELISA Kit
Reactivity	Human
Range	0.05-20ng/ml
Tested Applications	ELISA
Target	HSPA9
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	73,680
Uniprot ID	P38646
Solubility (25°C)	24ng/ml
Sensitivity	0.019ng/ml
Entrez	3313
Expiration Date	Please enquire.