

Product Datasheet

Human Probable RNA-binding Protein 24 ELISA Kit (orb1210080)

Description	Human Probable RNA-binding Protein 24, RBM24 ELISA Kit
Reactivity	Human
Range	2-700ng/L
Tested Applications	ELISA
Target	RBM24
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	60,023
Uniprot ID	Q8TBYO
Solubility (25°C)	800ng/L
Sensitivity	1.14ng/L
Entrez	166863
Expiration Date	Please enquire.