

Product Datasheet

Human Cerebral Dopamine Neurotrophic Factor ELISA Kit (orb1210092)

Description	Human Cerebral Dopamine Neurotrophic Factor, CDNF ELISA Kit
Reactivity	Human
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	CDNF
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	20,964
Uniprot ID	Q49AH0
Solubility (25°C)	48ng/ml
Sensitivity	0.052ng/ml
Entrez	441549
Expiration Date	Please enquire.