

## Product Datasheet

# Human Centromere-associated Protein E ELISA Kit (orb1210094)

|                            |  |
|----------------------------|--|
| <b>Description</b>         | Human Centromere-associated Protein E, CENPE ELISA Kit |
| <b>Reactivity</b>          | Human  |
| <b>Range</b>               | 0.1-35ng/ml  |
| <b>Tested Applications</b> | ELISA  |
| <b>Target</b>              | CENPE  |
| <b>Storage</b>             | 2-8°C  |
| <b>Note</b>                | For research use only                                  |
| <b>Application notes</b>   | Sandwich   |
| <b>Sample Types</b>        | Serum, plasma, cell culture supernates                 |
| <b>MW</b>                  | 316,415  |
| <b>Uniprot ID</b>          | <a href="#">Q02224</a>                                 |
| <b>Solubility (25°C)</b>   | 40ng/ml  |
| <b>Sensitivity</b>         | 0.049ng/ml   |
| <b>Entrez</b>              | <a href="#">1062</a>                                   |
| <b>Expiration Date</b>     | Please enquire.  |