

## Product Datasheet

# Human Exosome Complex Component Rrp41 ELISA Kit (orb1210169)

<b>Description</b>	Human Exosome Complex Component Rrp41, EXOSC4 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	3-900ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	EXOSC4
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	26,383
<b>Uniprot ID</b>	<a href="#">Q9NPD3</a>
<b>Solubility (25°C)</b>	960ng/L
<b>Sensitivity</b>	1.15ng/L
<b>Entrez</b>	<a href="#">54512</a>
<b>Expiration Date</b>	Please enquire.