

Product Datasheet

Human Exosome Complex Component Rrp40 ELISA Kit (orb1210170)

Description	Human Exosome Complex Component Rrp40, EXOSC3 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	EXOSC3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,572
Uniprot ID	Q9NQT5
Solubility (25°C)	1600ng/L
Sensitivity	3.32ng/L
Entrez	51010
Expiration Date	Please enquire.