

Product Datasheet

Human Exosome Complex Component Rrp40 ELISA Kit (orb1210170)

Description	Human Exosome Complex Component Rrp40, EXOSC3 ELISA Kit
Reactivity	Human
Range	7-1500ng/L
Tested Applications	ELISA
Target	EXOSC3
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	29,572
Uniprot ID	Q9NQ75
Solubility (25°C)	1600ng/L
Sensitivity	3.32ng/L
Entrez	51010
Expiration Date	Please enquire.

Biorbyt Ltd.

7 Signet Court, Swann's Road,
Cambridge, CB5 8LA, United Kingdom
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+44 \(0\) 1223 859-353](tel:+44(0)1223859353) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)

Biorbyt LLC.

68 TW Alexander Drive,
Durham, NC, 27713, United States
Email: info@biorbyt.com, support@biorbyt.com
Phone: [+1 \(415\) 906-5211](tel:+1(415)9065211) | Fax: [+1 \(415\) 651-8558](tel:+1(415)6518558)