

## Product Datasheet

# Human Exosome Complex Component Rrp4 ELISA Kit (orb1210171)

<b>Description</b>	Human Exosome Complex Component Rrp4, EXOSC2 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	5-1000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	EXOSC2
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	32,789
<b>Uniprot ID</b>	<a href="#">Q13868</a>
<b>Solubility (25°C)</b>	1280ng/L
<b>Sensitivity</b>	2.23ng/L
<b>Entrez</b>	<a href="#">23404</a>
<b>Expiration Date</b>	Please enquire.