

Product Datasheet

Human Zymogen Granule Protein 16 Homolog B ELISA Kit (orb1210207)

| | |
|----------------------------|---|
| Description | Human Zymogen Granule Protein 16 Homolog B, ZG16B ELISA Kit |
| Reactivity | Human |
| Range | 50-10000ng/L |
| Tested Applications | ELISA |
| Target | ZG16B |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 22,739 |
| Uniprot ID | Q96DA0 |
| Solubility (25°C) | 12800ng/L |
| Sensitivity | 25.39ng/L |
| Entrez | 124220 |
| Expiration Date | Please enquire. |