

Product Datasheet

Mouse C-C Motif Chemokine 20 ELISA Kit (orb1210209)

Description	Mouse C-C Motif Chemokine 20, CCL20 ELISA Kit
Reactivity	Mouse
Range	2-600ng/L
Tested Applications	ELISA
Target	CCL20
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	10,826
Uniprot ID	O89093
Solubility (25°C)	640ng/L
Sensitivity	1.25ng/L
Entrez	20297
Expiration Date	Please enquire.