

Product Datasheet

Rat Serine ELISA Kit (orb1210219)

Description	Rat Serine, Threonine-protein Kinase Tao2, TAOK2 ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	TAOK2
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	138,751
Uniprot ID	Q9JLS3
Solubility (25°C)	3200ng/L
Sensitivity	7.33ng/L
Entrez	64666
Expiration Date	Please enquire.