

Product Datasheet

Rat Erythroferrone ELISA Kit (orb1210368)

| | |
|----------------------------|--|
| Description | Rat Erythroferrone, ERFE ELISA Kit |
| Reactivity | Rat |
| Range | 20-3800ng/L |
| Tested Applications | ELISA |
| Target | ERFE |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 36,369 |
| Uniprot ID | D4AB34 |
| Solubility (25°C) | 4000ng/L |
| Sensitivity | 10.49ng/L |
| Entrez | 681056 |
| Expiration Date | Please enquire. |