

Product Datasheet

Human Ephrin Type A Receptor 10 ELISA Kit (orb1210371)

Description	Human Ephrin Type A Receptor 10, EPHA10 ELISA Kit
Reactivity	Human
Range	10-2000ng/L
Tested Applications	ELISA
Target	EPHA10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	109,716
Uniprot ID	Q5JZY3
Solubility (25°C)	2400ng/L
Sensitivity	4.47ng/L
Entrez	284656
Expiration Date	Please enquire.