

Product Datasheet

Mouse Repetin ELISA Kit (orb1210388)

Description	Mouse Repetin, RPTN ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	RPTN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	128,619
Uniprot ID	P97347
Solubility (25°C)	48ng/mL
Sensitivity	0.055ng/mL
Entrez	20129
Expiration Date	Please enquire.