

Product Datasheet

Mouse Protein Wnt-11 ELISA Kit (orb1210414)

Description	Mouse Protein Wnt-11, WNT11 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	WNT11
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,135
Uniprot ID	P48615
Solubility (25°C)	48ng/mL
Sensitivity	0.054ng/mL
Entrez	22411
Expiration Date	Please enquire.