

Product Datasheet

Mouse Protein Disulfide-isomerase A4 ELISA Kit (orb1210428)

Description	Mouse Protein Disulfide-isomerase A4, PDIA4 ELISA Kit
Reactivity	Mouse
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	PDIA4
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	71,982
Uniprot ID	P08003
Solubility (25°C)	48ng/mL
Sensitivity	0.058ng/mL
Entrez	12304
Expiration Date	Please enquire.