

Product Datasheet

Mouse Prostaglandin-H2 D-isomerase ELISA Kit (orb1210432)

Description	Mouse Prostaglandin-H2 D-isomerase, PTGDS ELISA Kit
Reactivity	Mouse
Range	2-600ng/mL
Tested Applications	ELISA
Target	PTGDS
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,066
Uniprot ID	O09114
Solubility (25°C)	640ng/mL
Sensitivity	0.98ng/mL
Entrez	19215
Expiration Date	Please enquire.