

Product Datasheet

Mouse Prosaas ELISA Kit (orb1210433)

Description	Mouse Prosaas, PCSK1N ELISA Kit
Reactivity	Mouse
Range	10-2000ng/L
Tested Applications	ELISA
Target	PCSK1N
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,270
Uniprot ID	Q9QXV0
Solubility (25°C)	2400ng/L
Sensitivity	5.59ng/L
Entrez	30052
Expiration Date	Please enquire.