

## Product Datasheet

# Human Usp6 N-terminal Like Protein ELISA Kit (orb1210458)

<b>Description</b>	Human Usp6 N-terminal Like Protein, USP6NL ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	20-3500ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	USP6NL
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	94,104
<b>Uniprot ID</b>	<a href="#">Q92738</a>
<b>Solubility (25°C)</b>	4000ng/L
<b>Sensitivity</b>	9.35ng/L
<b>Entrez</b>	<a href="#">9712</a>
<b>Expiration Date</b>	Please enquire.