

## Product Datasheet

# Rat Transient Receptor Potential Cation Channel Subfamily V Member 4 ELISA Kit (orb1210522)

<b>Description</b>	Rat Transient Receptor Potential Cation Channel Subfamily V Member 4, TRPV4 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	15-3000ng/L
<b>Tested Applications</b>	ELISA
<b>Target</b>	TRPV4
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	98,010
<b>Uniprot ID</b>	<a href="#">Q9ERZ8</a>
<b>Solubility (25°C)</b>	3200ng/L
<b>Sensitivity</b>	7.19ng/L
<b>Entrez</b>	<a href="#">66026</a>
<b>Expiration Date</b>	Please enquire.