

Product Datasheet

Rat Transferrin Receptor Protein 1 ELISA Kit (orb1210530)

Description	Rat Transferrin Receptor Protein 1, CD71 ELISA Kit
Reactivity	Rat
Range	10-2000ng/mL
Tested Applications	ELISA
Target	TFR
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	85,876
Uniprot ID	Q99376
Solubility (25°C)	2400ng/mL
Sensitivity	5.08ng/mL
Entrez	64678
Expiration Date	Please enquire.