

Product Datasheet

Rat T-cell Surface Glycoprotein Cd3 Delta Chain ELISA Kit (orb1210546)

Description	Rat T-cell Surface Glycoprotein Cd3 Delta Chain, CD3D ELISA Kit
Reactivity	Rat
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	CD3D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	19,356
Uniprot ID	P19377
Solubility (25°C)	48ng/mL
Sensitivity	0.056ng/mL
Entrez	25710
Expiration Date	Please enquire.