

Product Datasheet

Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 8 ELISA Kit (orb1210608)

Description	Human Carcinoembryonic Antigen Related Cell Adhesion Molecule 8, CEACAM8 ELISA Kit
Reactivity	Human
Range	0.3-90ng/mL
Tested Applications	ELISA
Target	CEACAM8
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	43,320
Uniprot ID	A8MTB9
Solubility (25°C)	96ng/mL
Sensitivity	0.13ng/mL
Entrez	729767
Expiration Date	Please enquire.