

Product Datasheet

Mouse L-lactate Dehydrogenase B Chain ELISA Kit (orb1210635)

Description	Mouse L-lactate Dehydrogenase B Chain, LDHB ELISA Kit
Reactivity	Mouse
Range	0.5-300mU/mL
Tested Applications	ELISA
Target	LDHB
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	36,572
Uniprot ID	P16125
Solubility (25°C)	320mU/mL
Sensitivity	0.25mU/mL
Entrez	16832
Expiration Date	Please enquire.