

Product Datasheet

Rat Proenkephalin-B ELISA Kit (orb1210639)

Description	Rat Proenkephalin-B, PDYN ELISA Kit
Reactivity	Rat
Range	10-2000ng/L
Tested Applications	ELISA
Target	PDYN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	28,079
Uniprot ID	P06300
Solubility (25°C)	2400ng/L
Sensitivity	5.27ng/L
Entrez	29190
Expiration Date	Please enquire.