

Product Datasheet

Bovine Decorin ELISA Kit (orb1210651)

Description	Bovine Decorin, DCN ELISA Kit
Reactivity	Bovine
Range	0.6-180ng/mL
Tested Applications	ELISA
Target	DCN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,879
Uniprot ID	P21793
Solubility (25°C)	192ng/mL
Sensitivity	0.31ng/mL
Entrez	280760
Expiration Date	Please enquire.