

Product Datasheet

Human Interleukin 24 ELISA Kit (orb1210678)

Description	Human Interleukin 24, IL-24 ELISA Kit
Reactivity	Human
Range	5-1500ng/L
Tested Applications	ELISA
Target	IL-24
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,825
Uniprot ID	Q13007
Solubility (25°C)	1600ng/L
Sensitivity	2.56ng/L
Entrez	11009
Expiration Date	Please enquire.