

Product Datasheet

Human Interleukin 17D ELISA Kit (orb1210679)

Description	Human Interleukin 17D, IL-17D ELISA Kit
Reactivity	Human
Range	0.2-60ng/L
Tested Applications	ELISA
Target	IL-17D
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	21,893
Uniprot ID	Q8TAD2
Solubility (25°C)	80ng/L
Sensitivity	0.11ng/L
Entrez	53342
Expiration Date	Please enquire.