

Product Datasheet

Mouse IgG Receptor Fc γ n Large Subunit P51 ELISA Kit (orb1210689)

Description	Mouse IgG Receptor Fc γ n Large Subunit P51, FCGR1 ELISA Kit
Reactivity	Mouse
Range	15-3000ng/L
Tested Applications	ELISA
Target	FCGR1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	40,093
Uniprot ID	Q61559
Solubility (25°C)	3200ng/L
Sensitivity	7.05ng/L
Entrez	14132
Expiration Date	Please enquire.