

Product Datasheet

Rat Lactase-Phlorizin Hydrolase ELISA Kit (orb1210699)

Description	Rat Lactase-Phlorizin Hydrolase, LCT ELISA Kit
Reactivity	Rat
Range	0.05-20ng/mL
Tested Applications	ELISA
Target	LCT
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	217,268
Uniprot ID	Q02401
Solubility (25°C)	24ng/mL
Sensitivity	0.019ng/mL
Entrez	116569
Expiration Date	Please enquire.