

Product Datasheet

Mouse Ghrelin O-Acyltransferase ELISA Kit (orb1210714)

| | |
|----------------------------|---|
| Description | Mouse Ghrelin O-Acyltransferase, MBOAT4 ELISA Kit |
| Reactivity | Mouse |
| Range | 3-900ng/L |
| Tested Applications | ELISA |
| Target | MBOAT4 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 49,879 |
| Uniprot ID | POC7A3 |
| Solubility (25°C) | 960ng/L |
| Sensitivity | 1.49ng/L |
| Entrez | 234155 |
| Expiration Date | Please enquire. |