

Product Datasheet

Human Receptor-binding Cancer Antigen Expressed On Siso Cells Rcas1 ELISA Kit (orb1210746)

Description	Human Receptor-binding Cancer Antigen Expressed On Siso Cells Rcas1 ELISA Kit
Reactivity	Human
Range	0.5-100U/ml
Tested Applications	ELISA
Target	RCAS1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,377
Uniprot ID	O00559
Solubility (25°C)	128U/ml
Sensitivity	0.28U/ml
Entrez	9166
Expiration Date	Please enquire.