

Product Datasheet

Rat Cytochrome P450 11B2 ELISA Kit (orb1210788)

| | |
|----------------------------|---|
| Description | Rat Cytochrome P450 11B2, Mitochondrial, CYP11B2 ELISA Kit |
| Reactivity | Rat |
| Range | 20-6000ng/L |
| Tested Applications | ELISA |
| Target | CYP11B2 |
| Storage | 2-8°C |
| Note | For research use only |
| Application notes | Sandwich |
| Sample Types | Serum, plasma, cell culture supernates |
| MW | 58,241 |
| Uniprot ID | P30099 |
| Solubility (25°C) | 6400ng/L |
| Sensitivity | 9.73ng/L |
| Entrez | 24294 |
| Expiration Date | Please enquire. |