

Product Datasheet

Rat Biglycan ELISA Kit (orb1210825)

Description	Rat Biglycan, BGN ELISA Kit
Reactivity	Rat
Range	0.1-40ng/mL
Tested Applications	ELISA
Target	BGN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	41,706
Uniprot ID	P47853
Solubility (25°C)	48ng/mL
Sensitivity	0.043ng/mL
Entrez	25181
Expiration Date	Please enquire.