

Product Datasheet

Mouse B- and T-lymphocyte Attenuator ELISA Kit (orb1210842)

Description	Mouse B- and T-lymphocyte Attenuator, BTLA ELISA Kit
Reactivity	Mouse
Range	0.3-90ng/ml
Tested Applications	ELISA
Target	BTLA
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	34,337
Uniprot ID	Q7TSA3
Solubility (25°C)	96ng/ml
Sensitivity	0.12ng/ml
Entrez	208154
Expiration Date	Please enquire.