

Product Datasheet

Mouse Decorin ELISA Kit (orb1210919)

Description	Mouse Decorin, DCN ELISA Kit
Reactivity	Mouse
Range	0.05-30ng/ml
Tested Applications	ELISA
Target	DCN
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	39,809
Uniprot ID	P28654
Solubility (25°C)	32ng/ml
Sensitivity	0.023ng/ml
Entrez	13179
Expiration Date	Please enquire.