

Product Datasheet

Rat Peroxiredoxin 6 ELISA Kit (orb1210930)

Description	Rat Peroxiredoxin 6, PRDX6 ELISA Kit
Reactivity	Rat
Range	10.23-3000ng/L
Tested Applications	ELISA
Target	PRDX6
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	24,819
Uniprot ID	O35244
Solubility (25°C)	3200ng/L
Sensitivity	5.17ng/L
Entrez	94167
Expiration Date	Please enquire.