

Product Datasheet

Mouse G Protein Coupled Estrogen Receptor 1 ELISA Kit (orb1210934)

Description	Mouse G Protein Coupled Estrogen Receptor 1, GPER ELISA Kit
Reactivity	Mouse
Range	0.2-70ng/ml
Tested Applications	ELISA
Target	GPER
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	42,479
Uniprot ID	Q8BMP4
Solubility (25°C)	80ng/ml
Sensitivity	0.13ng/ml
Entrez	76854
Expiration Date	Please enquire.