

Product Datasheet

Rat Transient Receptor Potential Cation Channel Subfamily V ELISA Kit (orb1210948)

Description	Rat Transient Receptor Potential Cation Channel Subfamily V, Member 1, TRPV1 ELISA Kit
Reactivity	Rat
Range	15-3000ng/L
Tested Applications	ELISA
Target	TRPV1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	94,948
Uniprot ID	O35433
Solubility (25°C)	3200ng/L
Sensitivity	7.13ng/L
Entrez	83810
Expiration Date	Please enquire.