

Product Datasheet

Rat Prostaglandin E Synthase ELISA Kit (orb1210962)

Description	Rat Prostaglandin E Synthase, PTGES ELISA Kit
Reactivity	Rat
Range	10-200ng/L
Tested Applications	ELISA
Target	PTGES
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	17,273
Uniprot ID	Q9JHF3
Solubility (25°C)	2400ng/L
Sensitivity	5.52ng/L
Entrez	59103
Expiration Date	Please enquire.