

Product Datasheet

Human Grb10-interacting Gyf Protein 1 ELISA Kit (orb1210977)

Description	Human Grb10-interacting Gyf Protein 1, GIGYF1 ELISA Kit
Reactivity	Human
Range	15-3000ng/L
Tested Applications	ELISA
Target	GIGYF1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	114,601
Uniprot ID	O75420
Solubility (25°C)	3200ng/L
Sensitivity	6.69ng/L
Entrez	64599
Expiration Date	Please enquire.