

Product Datasheet

Rat Microtubule Associated Protein Tau ELISA Kit (orb1211038)

Description	Rat Microtubule Associated Protein Tau, MAPT ELISA Kit
Reactivity	Rat
Range	5-1500ng/L
Tested Applications	ELISA
Target	MAPt
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	78,564
Uniprot ID	P19332
Solubility (25°C)	1920ng/L
Sensitivity	2.56ng/L
Entrez	29477
Expiration Date	Please enquire.