

## Product Datasheet

# Human A Disintegrin and Metalloproteinase with Thrombospondin Motifs 7 ELISA Kit (orb1211052)

<b>Description</b>	Human A Disintegrin and Metalloproteinase with Thrombospondin Motifs 7, ADAMTS7 ELISA Kit
<b>Reactivity</b>	Human
<b>Range</b>	1-400ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	ADAMTS7
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	184,095
<b>Uniprot ID</b>	<a href="#">Q9UKP4</a>
<b>Solubility (25°C)</b>	480ng/ml
<b>Sensitivity</b>	0.65ng/ml
<b>Entrez</b>	<a href="#">11173</a>
<b>Expiration Date</b>	Please enquire.