

## Product Datasheet

### Rat Kallikrein 6 ELISA Kit (orb1211053)

<b>Description</b>	Rat Kallikrein 6, KLK6 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.5-100ng/mL
<b>Tested Applications</b>	ELISA
<b>Target</b>	KLK6
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	29,013
<b>Uniprot ID</b>	<a href="#">P36374</a>
<b>Solubility (25°C)</b>	128ng/mL
<b>Sensitivity</b>	0.214ng/mL
<b>Entrez</b>	<a href="#">29245</a>
<b>Expiration Date</b>	Please enquire.