

Product Datasheet

Mouse Sclerostin Domain-containing Protein 1 ELISA Kit (orb1211056)

Description	Mouse Sclerostin Domain-containing Protein 1, SOSTDC1 ELISA Kit
Reactivity	Mouse
Range	3-900ng/L
Tested Applications	ELISA
Target	SOSTDC1
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	23,174
Uniprot ID	Q9CQN4
Solubility (25°C)	960ng/L
Sensitivity	1.53ng/L
Entrez	66042
Expiration Date	Please enquire.