

## Product Datasheet

# Rat A Disintegrin and Metalloproteinase with Thrombospondin Motifs 1 ELISA Kit (orb1211106)

<b>Description</b>	Rat A Disintegrin and Metalloproteinase with Thrombospondin Motifs 1, ADAMTS1 ELISA Kit
<b>Reactivity</b>	Rat
<b>Range</b>	0.05-15ng/ml
<b>Tested Applications</b>	ELISA
<b>Target</b>	ADAMTS1
<b>Storage</b>	2-8°C
<b>Note</b>	For research use only
<b>Application notes</b>	Sandwich
<b>Sample Types</b>	Serum, plasma, cell culture supernates
<b>MW</b>	105,706
<b>Uniprot ID</b>	<a href="#">Q9WUQ1</a>
<b>Solubility (25°C)</b>	16ng/ml
<b>Sensitivity</b>	0.027ng/ml
<b>Entrez</b>	<a href="#">79252</a>
<b>Expiration Date</b>	Please enquire.