

Product Datasheet

Rat Myelin-oligodendrocyte Glycoprotein ELISA Kit (orb1211108)

Description	Rat Myelin-oligodendrocyte Glycoprotein, MOG ELISA Kit
Reactivity	Rat
Range	0.1-40ng/ml
Tested Applications	ELISA
Target	MOG
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	27,882
Uniprot ID	Q63345
Solubility (25°C)	48ng/ml
Sensitivity	0.043ng/ml
Entrez	24558
Expiration Date	Please enquire.