

Product Datasheet

Human Kallikrein 10 ELISA Kit (orb1211120)

Description	Human Kallikrein 10, KLK10 ELISA Kit
Reactivity	Human
Range	0.1-50ng/ml
Tested Applications	ELISA
Target	KLK10
Storage	2-8°C
Note	For research use only
Application notes	Sandwich
Sample Types	Serum, plasma, cell culture supernates
MW	30,170
Uniprot ID	043240
Solubility (25°C)	56ng/ml
Sensitivity	0.06ng/ml
Entrez	5655
Expiration Date	Please enquire.